

PLD2 Antibody



PACO54938

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May have a role in signal-induced cytoskeletal regulation and/or endocytosis.

Gene ID:

PLD2

Uniprot

O14939

Synonyms:

Phospholipase D2 (PLD 2) (hPLD2) (EC 3.1.4.4) (Choline phosphatase 2) (PLD1C) (Phosphatidylcholine-hydrolyzing phospholipase D2), PLD2


Immunogen:

Recombinant Human Phospholipase D2 protein (1-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO54938 at dilution of 1:100.