

# SOX11 Antibody, FITC conjugated



PACO55020

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcriptional factor involved in the embryonic neurogenesis. May also have a role in tissue modeling during development.

**Gene ID:**

SOX11

**Uniprot**

P35716

**Synonyms:**

Transcription factor SOX-11, SOX11

**Immunogen:**

Recombinant Human Transcription factor SOX-11 protein (241-341AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A