

TOB1 Antibody, FITC conjugated



PACO55028

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Anti-proliferative protein; the function is mediated by association with deadenylase subunits of the CCR4-NOT complex. Mediates CPEB3-accelerated mRNA deadenylation by binding to CPEB3 and recruiting CNOT7 which leads to target mRNA deadenylation and decay.

Gene ID:

TOB1

Uniprot

P50616

Synonyms:

Protein Tob1 (Transducer of erbB-2 1), TOB1, TOB TROB1

Immunogen:

Recombinant Human Protein Tob1 protein (42-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A