

EC12 Antibody, HRP conjugated



PACO55091

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Able to isomerize both 3-cis and 3-trans double bonds into the 2-trans form in a range of enoyl-CoA species. Has a preference for 3-trans substrates.

Gene ID:

EC12

Uniprot

O75521

Synonyms:

Enoyl-CoA delta isomerase 2, mitochondrial (EC 5.3.3.8) (DRS-1) (Delta(3), delta(2)-enoyl-CoA isomerase) (D3, D2-enoyl-CoA isomerase) (Diazepam-binding inhibitor-related protein 1) (DBI-related protein 1) (Dodecenoyl-CoA isomerase) (Hepatocellular carcinoma-associated antigen 88) (Peroxisomal 3,2-trans-enoyl-CoA isomerase) (pECI) (Renal carcinoma antigen NY-REN-1), EC12, DRS1 HCA88 Peci

Immunogen:

Recombinant Human Enoyl-CoA delta isomerase 2, mitochondrial protein (119-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A