

GAS6 Antibody, Biotin conjugated



PACO55133

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ligand for tyrosine-protein kinase receptors AXL, TYRO3 and MER whose signaling is implicated in cell growth and survival, cell adhesion and cell migration. GAS6/AXL signaling plays a role in various processes such as endothelial cell survival during acidification by preventing apoptosis, optimal cytokine signaling during human natural killer cell development, hepatic regeneration, gonadotropin-releasing hormone neuron survival and migration, platelet activation, or regulation of thrombotic responses.

Gene ID:

GAS6

Uniprot

Q14393

Synonyms:

Growth arrest-specific protein 6 (GAS-6) (AXL receptor tyrosine kinase ligand), GAS6, AXLLG

Immunogen:

Recombinant Human Growth arrest-specific protein 6 protein (70-217AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A