

SLC16A12 Antibody, Biotin conjugated



PACO55189

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Proton-linked monocarboxylate transporter that mediates creatine transport across the plasma membrane.

Gene ID:

SLC16A12

Uniprot

Q6ZSM3

Synonyms:

Monocarboxylate transporter 12 (MCT 12) (Creatine transporter 2) (CRT2) (Solute carrier family 16 member 12), SLC16A12, MCT12

Immunogen:

Recombinant Human Monocarboxylate transporter 12 protein (431-486AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A