

SH3TC2 Antibody



PACO55214

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

SH3TC2

Uniprot

Q8TF17

Synonyms:

SH3 domain and tetratricopeptide repeat-containing protein 2, SH3TC2, KIAA1985


Immunogen:

Recombinant Human SH3 domain and tetratricopeptide repeat-containing protein 2 protein (188-427AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO55214 at dilution of 1:100.



Immunofluorescent analysis of PC-3 cells using PACO55214 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).