

PANX1 Antibody, Biotin conjugated



PACO55261

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Structural component of the gap junctions and the hemichannels. May play a role as a Ca(2+)-leak channel to regulate ER Ca(2+) homeostasis.

Gene ID:

PANX1

Uniprot

Q96RD7

Synonyms:

Pannexin-1, PANX1, MRS1

Immunogen:

Recombinant Human Pannexin-1 protein (288-426AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A