

FOXD3 Antibody, FITC conjugated



PACO55344

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds to the consensus sequence 5'-A[AT]T[AG]TTTGTTT-3' and acts as a transcriptional repressor. Also acts as a transcriptional activator. Promotes development of neural crest cells from neural tube progenitors. Restricts neural progenitor cells to the neural crest lineage while suppressing interneuron differentiation. Required for maintenance of pluripotent cells in the pre-implantation and peri-implantation stages of embryogenesis.

Gene ID:

FOXD3

Uniprot

Q9UJU5

Synonyms:

Forkhead box protein D3 (HNF3/FH transcription factor genesis), FOXD3, HFH2

Immunogen:

Recombinant Human Forkhead box protein D3 protein (1-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A