

ADH1A Antibody, HRP conjugated



PACO55383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, nucleoplasm, plasma membrane, alcohol dehydrogenase (NAD) activity, alcohol dehydrogenase activity, zinc-dependent, retinol dehydrogenase activity, alcohol metabolic process, drug metabolic process, ethanol oxidation.

Gene ID:

ADH1A

Uniprot

P07327

Synonyms:

Alcohol dehydrogenase 1A (EC 1.1.1.1) (Alcohol dehydrogenase subunit alpha), ADH1A, ADH1

Immunogen:

Recombinant Human Alcohol dehydrogenase 1A protein (104-253AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A