

BCL3 Antibody, Biotin conjugated



PACO55409

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Contributes to the regulation of transcriptional activation of NF-kappa-B target genes. In the cytoplasm, inhibits the nuclear translocation of the NF-kappa-B p50 subunit. In the nucleus, acts as transcriptional activator that promotes transcription of NF-kappa-B target genes. Contributes to the regulation of cell proliferation.

Gene ID:

BCL3

Uniprot

P20749

Synonyms:

B-cell lymphoma 3 protein (BCL-3) (Proto-oncogene BCL3), BCL3, BCL4 D19S37

Immunogen:

Recombinant Human B-cell lymphoma 3 protein (362-454AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A