

CREM Antibody, Biotin conjugated



PACO55441

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional regulator that binds the cAMP response element (CRE), a sequence present in many viral and cellular promoters. Isoforms are either transcriptional activators or repressors. Plays a role in spermatogenesis and is involved in spermatid maturation.

Gene ID:

CREM

Uniprot

Q03060

Synonyms:

cAMP-responsive element modulator (Inducible cAMP early repressor) (ICER), CREM

Immunogen:

Recombinant Human cAMP-responsive element modulator protein (145-321AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A