

CTNNA2 Antibody, HRP conjugated



PACO55443

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as a linker between cadherin adhesion receptors and the cytoskeleton to regulate cell-cell adhesion and differentiation in the nervous system. Regulates morphological plasticity of synapses and cerebellar and hippocampal lamination during development. Functions in the control of startle modulation.

Gene ID:

CTNNA2

Uniprot

P26232

Synonyms:

Catenin alpha-2 (alpha N-catenin) (Alpha-catenin-related protein), CTNNA2, CAPR

Immunogen:

Recombinant Human Catenin alpha-2 protein (826-939AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A