

EOMES Antibody, FITC conjugated



PACO55476

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as a transcriptional activator playing a crucial role during development. Functions in trophoblast differentiation and later in gastrulation, regulating both mesoderm delamination and endoderm specification. Plays a role in brain development being required for the specification and the proliferation of the intermediate progenitor cells and their progeny in the cerebral cortex. Also involved in the differentiation of CD8+ T-cells during immune response regulating the expression of lytic effector genes.

Gene ID:

EOMES

Uniprot

O95936

Synonyms:

Eomesodermin homolog (T-box brain protein 2) (T-brain-2) (TBR-2), EOMES, TBR2

Immunogen:

Recombinant Human Eomesodermin homolog protein (562-676AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A