

EXOSC10 Antibody, FITC conjugated



PACO55484

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

EXOSC10

Uniprot

Q01780

Synonyms:

Exosome component 10 (EC 3.1.13) (Autoantigen PM/Scl 2) (P100 polymyositis-scleroderma overlap syndrome-associated autoantigen) (Polymyositis/scleroderma autoantigen 100 kDa) (PM/Scl-100) (Polymyositis/scleroderma autoantigen 2), EXOSC10, PMSCL PMSCL2 RRP6

Immunogen:

Recombinant Human Exosome component 10 protein (587-762AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A