

GGCX Antibody, HRP conjugated



PACO55523

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates the vitamin K-dependent carboxylation of glutamate residues to calcium-binding gamma-carboxyglutamate (Gla) residues with the concomitant conversion of the reduced hydroquinone form of vitamin K to vitamin K epoxide.

Gene ID:

GGCX

Uniprot

P38435

Synonyms:

Vitamin K-dependent gamma-carboxylase (EC 4.1.1.90) (Gamma-glutamyl carboxylase) (Peptidyl-glutamate 4-carboxylase) (Vitamin K gamma glutamyl carboxylase), GG CX, GC

Immunogen:

Recombinant Human Vitamin K-dependent gamma-carboxylase protein (383-526AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A