

GPR12 Antibody, HRP conjugated



PACO55527

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Promotes neurite outgrowth and blocks myelin inhibition in neurons. Receptor with constitutive G(s) signaling activity that stimulates cyclic AMP production.

Gene ID:

GPR12

Uniprot

P47775

Synonyms:

G-protein coupled receptor 12, GPR12

Immunogen:

Recombinant Human G-protein coupled receptor 12 protein (1-48AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A