

# MSN Antibody, FITC conjugated



PACO55568

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probably involved in connections of major cytoskeletal structures to the plasma membrane. May inhibit herpes simplex virus 1 infection at an early stage. Plays a role in regulating the proliferation, migration, and adhesion of human lymphoid cells and participates in immunologic synapse formation.

**Gene ID:**

MSN

**Uniprot**

P26038

**Synonyms:**

Moesin (Membrane-organizing extension spike protein), MSN

**Immunogen:**

Recombinant Human Moesin protein (466-572AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A