

MSN Antibody, Biotin conjugated



PACO55569

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probably involved in connections of major cytoskeletal structures to the plasma membrane. May inhibit herpes simplex virus 1 infection at an early stage. Plays a role in regulating the proliferation, migration, and adhesion of human lymphoid cells and participates in immunologic synapse formation.

Gene ID:

MSN

Uniprot

P26038

Synonyms:

Moesin (Membrane-organizing extension spike protein), MSN

Immunogen:

Recombinant Human Moesin protein (466-572AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A