

KCNAB1 Antibody, HRP conjugated



PACO55575

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KCNAB1

Uniprot

Q14722

Synonyms:

Voltage-gated potassium channel subunit beta-1 (EC 1.1.1) (K(+)) channel subunit beta-1) (Kv-beta-1), KCNAB1, KCNA1B

Immunogen:

Recombinant Human Voltage-gated potassium channel subunit beta-1 protein (287-401AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A