

MAP3K4 Antibody, Biotin conjugated



PACO55613

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of a protein kinase signal transduction cascade. Activates the CSBP2, P38 and JNK MAPK pathways, but not the ERK pathway. Specifically phosphorylates and activates MAP2K4 and MAP2K6.

Gene ID:

MAP3K4

Uniprot

Q9Y6R4

Synonyms:

Mitogen-activated protein kinase kinase kinase 4 (EC 2.7.11.25) (MAP three kinase 1) (MAPK/ERK kinase kinase 4) (MEK kinase 4) (MEKK 4), MAP3K4, KIAA0213 MAPKKK4 MEKK4 MTK1

Immunogen:

Recombinant Human Mitogen-activated protein kinase kinase kinase 4 protein (53-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A