

# MAP3K6 Antibody, FITC conjugated



PACO55616

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Component of a protein kinase signal transduction cascade. Activates the JNK, but not ERK or p38 kinase pathways.

**Gene ID:**

MAP3K6

**Uniprot**

O95382

**Synonyms:**

Mitogen-activated protein kinase kinase kinase 6 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 2), MAP3K6, ASK2 MAPKKK6 MEKK6

**Immunogen:**

Recombinant Human Mitogen-activated protein kinase kinase kinase 6 protein (905-990AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A