

# MC1R Antibody, HRP conjugated



PACO55627

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for MSH (alpha, beta and gamma) and ACTH. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

**Gene ID:**

MC1R

**Uniprot**

Q01726

**Synonyms:**

Melanocyte-stimulating hormone receptor (MSH-R) (Melanocortin receptor 1) (MC1-R), MC1R, MSHR

**Immunogen:**

Recombinant Human Melanocyte-stimulating hormone receptor protein (1-37AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A