

MCHR2 Antibody, FITC conjugated



PACO55632

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for melanin-concentrating hormone, coupled to G proteins that activate phosphoinositide hydrolysis.

Gene ID:

MCHR2

Uniprot

Q969V1

Synonyms:

Melanin-concentrating hormone receptor 2 (MCH receptor 2) (MCH-R2) (MCHR-2) (G-protein coupled receptor 145) (GPRv17) (MCH-2R) (MCH2) (MCH2R), MCHR2, GPR145 SLT

Immunogen:

Recombinant Human Melanin-concentrating hormone receptor 2 protein (1-39AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A