

MGST1 Antibody, HRP conjugated



PACO55643

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Has a wide substrate specificity.

Gene ID:

MGST1

Uniprot

P10620

Synonyms:

Microsomal glutathione S-transferase 1 (Microsomal GST-1) (EC 2.5.1.18) (Microsomal GST-I), MGST1, GST12 MGST

Immunogen:

Recombinant Human Microsomal glutathione S-transferase 1 protein (63-96AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A