

MYH7 Antibody, HRP conjugated



PACO55667

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Muscle contraction.

Gene ID:

MYH7

Uniprot

P12883

Synonyms:

Myosin-7 (Myosin heavy chain 7) (Myosin heavy chain slow isoform) (MyHC-slow) (Myosin heavy chain, cardiac muscle beta isoform) (MyHC-beta), MYH7, MYHCB

Immunogen:

Recombinant Human Myosin-7 protein (1802-1915AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A