

NOVA1 Antibody, HRP conjugated



PACO55695

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May regulate RNA splicing or metabolism in a specific subset of developing neurons.

Gene ID:

NOVA1

Uniprot

P51513

Synonyms:

RNA-binding protein Nova-1 (Neuro-oncological ventral antigen 1) (Onconeural ventral antigen 1) (Paraneoplastic Ri antigen) (Ventral neuron-specific protein 1), NOVA1

Immunogen:

Recombinant Human RNA-binding protein Nova-1 protein (6-174AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A