

# PEX10 Antibody



PACO55722

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

Somewhat implicated in the biogenesis of peroxisomes.

**Gene ID:**

PEX10

**Uniprot**

O60683

**Synonyms:**

Peroxisome biogenesis factor 10 (Peroxin-10) (Peroxisomal biogenesis factor 10) (Peroxisome assembly protein 10) (RING finger protein 69), PEX10, RNF69

**Immunogen:**


Recombinant Human Peroxisome biogenesis factor 10 protein (63-157AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Rat brain tissue, Mouse brain tissue. All lanes: PEX10 antibody at 2.7 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 38, 40 kDa. Observed band size: 38 kDa.