

# PEX10 Antibody, Biotin conjugated



PACO55725

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Somewhat implicated in the biogenesis of peroxisomes.

**Gene ID:**

PEX10

**Uniprot**

O60683

**Synonyms:**

Peroxisome biogenesis factor 10 (Peroxin-10) (Peroxisomal biogenesis factor 10) (Peroxisome assembly protein 10) (RING finger protein 69), PEX10, RNF69

**Immunogen:**

Recombinant Human Peroxisome biogenesis factor 10 protein (63-157AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A