

# POLD1 Antibody, FITC conjugated



PACO55740

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POLD1

**Uniprot**

P28340

**Synonyms:**

DNA polymerase delta catalytic subunit (EC 2.7.7.7) (EC 3.1.11) (DNA polymerase subunit delta p125), POLD1, POLD

**Immunogen:**

Recombinant Human DNA polymerase delta catalytic subunit protein (1-123AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A