

POLD1 Antibody, Biotin conjugated



PACO55741

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

POLD1

Uniprot

P28340

Synonyms:

DNA polymerase delta catalytic subunit (EC 2.7.7.7) (EC 3.1.11) (DNA polymerase subunit delta p125), POLD1, POLD

Immunogen:

Recombinant Human DNA polymerase delta catalytic subunit protein (1-123AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A