

RPS7 Antibody, HRP conjugated



PACO55763

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for rRNA maturation.

Gene ID:

RPS7

Uniprot

P62081

Synonyms:

40S ribosomal protein S7 (Small ribosomal subunit protein eS7), RPS7

Immunogen:

Recombinant Human 40S ribosomal protein S7 protein (1-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A