

TCP1 Antibody, HRP conjugated



PACO55803

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Molecular chaperone; assists the folding of proteins upon ATP hydrolysis. As part of the BBS/CCT complex may play a role in the assembly of BBSome, a complex involved in ciliogenesis regulating transports vesicles to the cilia. Known to play a role, in vitro, in the folding of actin and tubulin.

Gene ID:

TCP1

Uniprot

P17987

Synonyms:

T-complex protein 1 subunit alpha (TCP-1-alpha) (CCT-alpha), TCP1, CCT1 CCTA

Immunogen:

Recombinant Human T-complex protein 1 subunit α protein (406-506AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A