

# PNMT Antibody, HRP conjugated



PACO55843

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Converts noradrenaline to adrenaline.

**Gene ID:**

PNMT

**Uniprot**

P11086

**Synonyms:**

Phenylethanolamine N-methyltransferase (PNMTase) (EC 2.1.1.28) (Noradrenaline N-methyltransferase), PNMT, PENT

**Immunogen:**

Recombinant Human Phenylethanolamine N-methyltransferase protein (1-176AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A