

EPS8 Antibody, FITC conjugated



PACO55892

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EPS8

Uniprot

Q12929

Synonyms:

Epidermal growth factor receptor kinase substrate 8, EPS8

Immunogen:

Recombinant Human Epidermal growth factor receptor kinase substrate 8 protein (572-670AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A