

NFATC3 Antibody, Biotin conjugated



PACO55897

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a regulator of transcriptional activation. Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2.

Gene ID:

NFATC3

Uniprot

Q12968

Synonyms:

Nuclear factor of activated T-cells, cytoplasmic 3 (NF-ATc3) (NFATc3) (NFATx) (T-cell transcription factor NFAT4) (NF-AT4), NFATC3, NFAT4

Immunogen:

Recombinant Human Nuclear factor of activated T-cells, cytoplasmic 3 protein (887-1070AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A