

LAMB3 Antibody, FITC conjugated



PACO55900

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

Gene ID:

LAMB3

Uniprot

Q13751

Synonyms:

Laminin subunit beta-3 (Epiligrin subunit bata) (Kalinin B1 chain) (Kalinin subunit beta) (Laminin B1k chain) (Laminin-5 subunit beta) (Nicein subunit beta), LAMB3, LAMNB1

Immunogen:

Recombinant Human Laminin subunit beta-3 protein (256-397AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A