

MECOM Antibody, HRP conjugated



PACO55907

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Functions as a transcriptional regulator binding to DNA sequences in the promoter region of target genes and regulating positively or negatively their expression. Oncogene which plays a role in development, cell proliferation and differentiation. May also play a role in apoptosis through regulation of the JNK and TGF-beta signaling. Involved in hematopoiesis.

Gene ID:

MECOM

Uniprot

Q03112

Synonyms:

MDS1 and EVI1 complex locus protein (Ecotropic virus integration site 1 protein homolog) (EVI-1) (Myelodysplasia syndrome 1 protein) (Myelodysplasia syndrome-associated protein 1), MECOM, EVI1

Immunogen:

Recombinant Human Histone-lysine N-methyltransferase MECOM protein (534-632AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A