

UBAP2L Antibody, HRP conjugated



PACO55931

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in the activity of long-term repopulating hematopoietic stem cells (LT-HSCs).

Gene ID:

UBAP2L

Uniprot

Q14157

Synonyms:

Ubiquitin-associated protein 2-like (Protein NICE-4), UBAP2L, KIAA0144 NICE4

Immunogen:

Recombinant Human Ubiquitin-associated protein 2-like protein (104-252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A