

TRDMT1 Antibody, Biotin conjugated



PACO55949

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Specifically methylates cytosine 38 in the anticodon loop of tRNA(Asp).

Gene ID:

TRDMT1

Uniprot

O14717

Synonyms:

tRNA (cytosine(38)-C(5))-methyltransferase (EC 2.1.1.204) (DNA (cytosine-5)-methyltransferase-like protein 2) (Dnmt2) (DNA methyltransferase homolog HsaIIIP) (DNA MTase homolog HsaIIIP) (M. HsaIIIP) (PuMet), TRDMT1, DNMT2

Immunogen:

Recombinant Human tRNA (cytosine(38)-C(5))-methyltransferase protein (170-269AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A