

# DPP9 Antibody, HRP conjugated



PACO55991

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Dipeptidyl peptidase that cleaves off N-terminal dipeptides from proteins having a Pro or Ala residue at position 2.

**Gene ID:**

DPP9

**Uniprot**

Q86T12

**Synonyms:**

Dipeptidyl peptidase 9 (DP9) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 2) (DPRP-2) (Dipeptidyl peptidase IX) (DPP IX) (Dipeptidyl peptidase-like protein 9) (DPLP9), DPP9, DPRP2

**Immunogen:**

Recombinant Human Dipeptidyl peptidase 9 protein (289-437AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A