

# PNPLA6 Antibody, Biotin conjugated



PACO56013

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Phospholipase B that deacylates intracellular phosphatidylcholine (PtdCho), generating glycerophosphocholine (GroPtdCho). This deacylation occurs at both sn-2 and sn-1 positions of PtdCho. Its specific chemical modification by certain organophosphorus (OP) compounds leads to distal axonopathy.

**Gene ID:**

PNPLA6

**Uniprot**

Q8IY17

**Synonyms:**

Neuropathy target esterase (EC 3.1.1.5) (Patatin-like phospholipase domain-containing protein 6), PNPLA6, NTE

**Immunogen:**

Recombinant Human Neuropathy target esterase protein (345-493AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A