

GPBAR1 Antibody, FITC conjugated



PACO56024

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for bile acid, Bile acid, binding induces its internalization, activation of extracellular signal-regulated kinase and intracellular cAMP production. May be involved in the suppression of macrophage functions by bile acid, .

Gene ID:

GPBAR1

Uniprot

Q8TDU6

Synonyms:

G-protein coupled bile acid, receptor 1 (G-protein coupled receptor GPCR19) (hGPCR19) (Membrane-type receptor for bile acid,) (M-BAR) (hBG37) (BG37), GPBAR1, TGR5

Immunogen:

Recombinant Human G-protein coupled bile acid, receptor 1 protein (283-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A