

# GPBAR1 Antibody, Biotin conjugated



PACO56025

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for bile acid, Bile acid, binding induces its internalization, activation of extracellular signal-regulated kinase and intracellular cAMP production. May be involved in the suppression of macrophage functions by bile acid, .

**Gene ID:**

GPBAR1

**Uniprot**

Q8TDU6

**Synonyms:**

G-protein coupled bile acid, receptor 1 (G-protein coupled receptor GPCR19) (hGPCR19) (Membrane-type receptor for bile acid, ) (M-BAR) (hBG37) (BG37), GPBAR1, TGR5

**Immunogen:**

Recombinant Human G-protein coupled bile acid, receptor 1 protein (283-330AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A