

NUP93 Antibody, Biotin conjugated



PACO56041

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. May anchor nucleoporins, but not NUP153 and TPR, to the NPC. During renal development, regulates podocyte migration and proliferation through SMAD4 signaling.

Gene ID:

NUP93

Uniprot

Q8N1F7

Synonyms:

Nuclear pore complex protein Nup93 (93 kDa nucleoporin) (Nucleoporin Nup93), NUP93, KIAA0095

Immunogen:

Recombinant Human Nuclear pore complex protein Nup93 protein (14-180AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A