

CNNM3 Antibody, Biotin conjugated



PACO56061

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable metal transporter.

Gene ID:

CNNM3

Uniprot

Q8NE01

Synonyms:

Metal transporter CNNM3 (Ancient conserved domain-containing protein 3) (Cyclin-M3), CNNM3, ACDP3

Immunogen:

Recombinant Human Metal transporter CNNM3 protein (358-645AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A