

ARX Antibody, Biotin conjugated



PACO56073

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor required for normal brain development. May be important for maintenance of specific neuronal subtypes in the cerebral cortex and axonal guidance in the floor plate.

Gene ID:

ARX

Uniprot

Q96QS3

Synonyms:

Homeobox protein ARX (Aristaless-related homeobox), ARX

Immunogen:

Recombinant Human Homeobox protein ARX protein (381-460AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A