

PKP2 Antibody, HRP conjugated



PACO56119

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in junctional plaques.

Gene ID:

PKP2

Uniprot

Q99959

Synonyms:

Plakophilin-2, PKP2

Immunogen:

Recombinant Human Plakophilin-2 protein (61-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A