

# PKP2 Antibody, FITC conjugated



PACO56120

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in junctional plaques.

**Gene ID:**

PKP2

**Uniprot**

Q99959

**Synonyms:**

Plakophilin-2, PKP2

**Immunogen:**

Recombinant Human Plakophilin-2 protein (61-183AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A