

RNF19A Antibody, HRP conjugated



PACO56175

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

E3 ubiquitin-protein ligase which accepts ubiquitin from E2 ubiquitin-conjugating enzymes UBE2L3 and UBE2L6 in the form of a thioester and then directly transfers the ubiquitin to targeted substrates, such as SNCAIP or CASR. Specifically ubiquitinates pathogenic SOD1 variants, which leads to their proteasomal degradation and to neuronal protection.

Gene ID:

RNF19A

Uniprot

Q9NV58

Synonyms:

E3 ubiquitin-protein ligase RNF19A (EC 2.3.2.27) (Double ring-finger protein) (Dorfin) (RING finger protein 19A) (p38), RNF19A, RNF19

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase RNF19A protein (741-825AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A